


HOW TO CHOOSE YOUR **BTM** SYSTEM?



CAMERA	BTM-T5 BTM-TIC500 BTM-TIC600	Fever Detection, Over Temperature Alert and Face Recognition	Multiple Targets Simultaneously	Requires a Professional Installation	VMS and NVR	High	±0.3°C	Public Spaces Hospitals, Factories, Mining, Hotels, Clubs, Airports, Stations
ACCESS CONTROL TERMINAL	BTM-213H BTM-213W	Fever Detection, Over Temperature Alert, Face Recognition and Access Control	Individual Target	Requires a Professional Installation	Access Control System	Moderate	±0.3°C	Limited Public Access Private Offices, Control Rooms, Back Offices
STAND ALONE UNIT	BTM-CW180 BTM-OTC513	Fever Detection and Over Temperature Alert	Individual Target	DIY	No Integration	Low	±0.3°C	Controlled Public Access Restaurants, Nursing Homes, Schools
HANDHELD IMAGER	BTM-16K BTM-16H BTM-30R	Fever Detection and Over Temperature Alert	Individual Target	Not Required	No Integration	Lowest	±0.3°C	All Purpose

HOW TO CHOOSE YOUR BTM SYSTEM?


CAMERA



STAND ALONE UNIT










ACCESS CONTROL TERMINAL



HANDHELD IMAGER



HOW TO CHOOSE YOUR BTM SYSTEM?

CAMERA		ACCESS CONTROL TERMINAL		STAND ALONE UNIT		HANDHELD IMAGER	
PRODUCT		PRODUCT		PRODUCT		PRODUCT	
BTM-T5		BTM-213H		BTM-CW180		BTM-16H	
							
BTM-TIC500		BTM-213W		BTM-OTC513		BTM-16K	
							
BTM-TIC600						BTM-30R	
